

LABSOLUTE® SENSITIVE, NITRILE, POWDER FREE

Article numbers: 7.696 904 (XS), 7.696 900 (S), 7.696 901 (M), 7.696 902 (L),
7.696 903 (XL)

General information:

Type: disposable gloves / examination gloves
Description: All information is printed on the dispenser box.
Shape: ambidextrous, straight fingers
Material: Nitrile Butadiene Rubber (NBR)
Color: lavender blue (valid for HACCP Standard)
Inner: online chlorinated, powder free
Outer: no treatment
Cuff/Surface: rolled cuff, textured fingers
Length: Ø 240 mm
Thickness (double): Finger min. 0.14 mm, Palm min. 0.12 mm, Cuff min. 0.10 mm
Gloves width: median (±10 mm): XS <80 mm, S 80 mm, M 95 mm, L 110 mm, XL >110 mm
Force at break (EN 455-2): Ø 6,0 N (during shelf life)
Residual powder: ≤ 2 mg / gloves
Available sizes: XS (5 - 6), S (6 - 7), M (7 - 8), L (8 - 9), XL (9 - 10)

Storage:

Storage life: 3 years (expiry date is printed on the dispenser box).
Storage conditions: cool (max. 30 °C), dry, dust free, avoid UV radiation, contact with oxidizing agents, do not store close to photocopy equipment.
Miscellaneous: Copper ions can discolour the gloves.

Directives and Certification:

Directives: Regulation (EU) 2016/425, CE category III
Standards: EN 420:2003+A1:2009, ISO 374-1:2016+A1:2018 (type B), EN 374-2:2014 (Level 2 = AQL 1,5), EN 16523-1:2015+A1:2018, EN 374-4:2013, EN ISO 374-5:2016, EN455-1, -2, -3
Food regulation: Regulation (EU) No. 1935/2004

Warning:

Gloves contain the following allergenic substances: Dithiocarbamate, Zinc-Mercaptobenzothiazol
This product contains allergenic substances which can cause allergic reactions, including anaphylactic responses.