

LABSOLUTE® PROTECTIVE, NITRILE, POWDER FREE

Article numbers: 7.696 906 (S), 7.696 907 (M), 7.696 908 (L), 7.696 909 (XL)

General information:

Type: disposable gloves
Shape: ambidextrous, five straight fingers
Material: Nitrile Butadiene Rubber (NBR)
Color: blue (recommended for HACCP standard)
Inner: online chlorinated, powder free
Outer: no treatment
Cuff/Surface: rolled cuff, fully textured surface
Length: Ø 240 mm (EN 455-2)
Thickness (double): Finger min. 0.30 mm, Palm min. 0.22 mm, Cuff min. 0.20 mm
Force at break (EN 455-2): Ø 6,0 N (during shelf life)
Tensile strength/elongation: 14 MPa / 500 % (after aging: 14 MPa / 400 %)
Residual powder: ≤ 2 mg / gloves (EN 455-3)

Available sizes: S (6 - 7), M (7 - 8), L (8 - 9), XL (9 - 10)
Packing quantity: 100 pieces (X, M, L), 90 pieces (XL) per dispenser box

Storage:

Storage life: 3 years (expiry date is printed on the dispenser box)
Storage conditions: cool & dry (recommended: 10 - 30 °C at 60 - 80% humidity), dust free, avoid UV radiation, avoid contact with oxidizing agents and ozone

Directives and Certification:

Directives: (EC) 2016/425 Personal Protective Equipment, CE category III
Standards: EN 420:2003+A1:2009, ISO 374-1:2016+A1:2018 (Typ B), EN 374-2:2014 (Level 2 = AQL 1,5), EN 16523-1:2015+A1:2018, EN 374-4:2013, EN ISO 374-5:2016, ISO 2859-1:1999+A1:2011, EN 421:2010
Food regulation: Directive (EC) No. 1935/2004

Warning:

Gloves contain the following allergenic substances: Dithiocarbamate

This product contains allergenic substances which can cause allergic reactions, including anaphylactic responses.