

LABSOLUTE® DISPOSABLE GLOVES, LATEX, POWDER FREE

Article numbers: 7.696 911 (S), 7.696 912 (M), 7.696 913 (L), 7.696 914 (XL)

General information:

Type: disposable gloves / examination gloves
Description: All information is printed on the dispenser box.
Shape: ambidextrous, straight fingers
Material: natural rubber latex, low in allergens, powder free (protein content < 10 µg / g gloves)
Color: natural white
Inner: online chlorinated, powder free
Outer: no treatment
Cuff/Surface: rolled cuff, finger textured
Length: Ø 240 mm
Thickness (double): Finger min. 0.24 mm, Palm min. 0.20 mm, Cuff min. 0.12 mm
Force at break (EN 455-2): Ø 6.0 N (during shelf life)
Residual powder: ≤ 2 mg / gloves

Available sizes: S (6 - 7), M (7 - 8), L (8 - 9), XL (9 - 10)

Storage:

Storage life: 3 years (expiry date is printed on the dispenser box).
Storage conditions: cool (max. 30 °C), dry, dust free, avoid UV radiation, contact with oxidizing agents, do not store close to photocopy equipment.
Miscellaneous: Copper ions can discolour the gloves.

Directives and Certification:

Directives: (EC) 2016/425 PSA, Personal Protective Equipment, CE category III

Standards: EN 420:2003+A1:2009, ISO 374-1:2016/AMD 1:2018, EN 374-2:2014 (AQL value 1.5), EN 16523-1:2015+A1:2018, EN 374-4:2013, EN ISO 374-5:2016, ISO 2859-1:1999 AMD 1:2011
Manufactured according to EN 455.

Food regulation: 1935/2004

Warning:

Gloves contain the following allergenic substances: Dithiocarbamate, Zinc-Mercaptobenzothiazol

This product contains natural rubber latex and the mentioned allergenic substances which can cause allergic reactions, including anaphylactic responses.