

HAND PROTECTOR GRIP, HEATGUARD

The red LABSOLUTE[®] hand protector grip HeatGuard made from silicone rubber protects your hands best against hot and cold surfaces. You can store the hand protector grip by using one of the two eyelets.



Your advantages:

- Temperature resistant from about -50 up to +250 °C.
- Fits all hand sizes.
- Nobs on the HeatGuard make sure that there is a safe, slip resistant grip all the time.

Please send a message to pmlabsolute@thgeyer.de or ask your sales representative, if there are any questions about the chemical resistance of the powder funnels.

| Item no. | Description | Pieces per pack |
|-----------|--------------------------------|-----------------|
| 7 696 925 | Heat Protector Grip, Heatguard | 1 |