

Technical Datasheet

Article No. 7753

Water Bath Preservative (500x)

For the disinfection of different types of water baths, effectively prevents microbial growth (fungi and bacteria). Active ingredient: Bronopol.

Application:

Use 2 ml Waterbath-Preservative per liter of water. Every 1 – 2 weeks add 2 ml Water Bath Preservative per liter (degradation of the active ingredient). Every 4-6 weeks, the water should be completely renewed.

For Laboratory use.

Parameter	Specification
Appearance	clear solution
Sterility	complies
Storage	-20 °C