

Technical Datasheet

Article No. 800

Silver chloride
min. 74.8 % Ag

For laboratory use.

Parameter	Value
Appearance	White crystal solid
Aluminium (Al)	max. 0.001 %
Bismuth (Bi)	max. 0.001 %
Calcium (Ca)	max. 0.003 %
Cadmium (Cd)	max. 0.002 %
Copper (Cu)	max. 0.0015 %
Iron (Fe)	max. 0.001 %
Magnesium (Mg)	max. 0.003 %
Manganese (Mn)	max. 0.001 %
Nickel (Ni)	max. 0.001 %
Lead (Pb)	max. 0.002 %
Silicium (Si)	max. 0.0025 %
Tin (Sn)	max. 0.001 %
Titanium (Ti)	max. 0.001 %