

# Technical Datasheet

Article No. 8096

**Tryptone  
(Powder)**

For Laboratory use.

Parameter	Specification	Result
Amino Nitrogen (AN)	min. 3.9 %	
Total Nitrogen (TN)	min. 10 %	
Loss on drying	max. 6 %	
Ash	max. 15 %	
pH (2 %)	6.5–7.5	
Standard Plate Count	max. 5,000 col/g	
Molds and yeast	max. 100 col/g	
Coliforms	negative	
Salmonella	negative	