

# Technical Datasheet

---

**Article no.: 836**

**Hydrochloric acid p. a. (min. 37.0 %)**

HCl

For laboratory use.

Parameter	Value
Ammonium (NH <sub>4</sub> )	max. 0.0003 %
Arsenic (As)	max. 0.000005 %
Cadmium (Cd)	max. 0.000001 %
Copper (Cu)	max. 0.000005 %
Free chlorine (Cl <sub>2</sub> )	max. 0.00005 %
Iron (Fe)	max. 0.00005 %
Lead (Pb)	max. 0.000005 %
Residue on ignition	max. 0.0005 %
Sulphate (SO <sub>4</sub> )	max. 0.0001 %
Sulphite (SO <sub>3</sub> )	max. 0.0002 %
Zinc (Zn)	max. 0.0001 %

Version no. 02-2024-03-28