

Technical Datasheet

Article No. 8622

Di-Natriumhydrogenphosphat-Dihydrat p. a.
min. 99,5 % $\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$

For laboratory use.

Parameter	Value
Appearance	White, cristalline powder
Mol.-Weight	177.99 g/mol
Density (20 °C)	2.10 g/ml
Ammonium (NH_4)	max. 0.001 %
Arsenic (As)	max. 0.0001 %
Iron (Fe)	max. 0.0005 %
Heavy metals (as Pb)	max. 0.001 %
Chloride (Cl)	max. 0.001 %
Sulfate (SO_4)	max. 0.005 %
Loss on drying	20,0–20,4 %