

# Technical Datasheet

**Article No. 8664**

**Sodium phosphate dibasic (dihydrate) p. a.**  
**min. 99.5 % Na<sub>2</sub>HPO<sub>4</sub> x 2H<sub>2</sub>O**

**For Laboratory use.**

Parameter	Value
Assay	min. 99.5 %
pH (5.0 % water)	8.8–9.4
Chloride (Cl)	max. 0.001 %
Sulphate (SO <sub>4</sub> )	max. 0.005 %
Total nitrogen (N)	max. 0.001 %
Heavy metals (as Pb)	max. 0.001 %
Iron (Fe)	max. 0.001 %
Potassium (K)	max. 0.005 %
Mol.- Weight	177.99 g/mol