

TECHNICAL DATASHEET

Art.-No.: 9213

Yeast extract for microbiology

For laboratory use, ACS, ISO, Reag. Ph. Eur.

Parameter	Specification
Appearance	Light brown powder
Density	0.5 g/ml
Assay (ex N on dried sample)	> 60.0 %
Phosphorus (P)	< 2.5 %
Coagulation protein	Not detected
Heavy metals (as Pb)	< 0.005 %
Loss on drying at 105 °C	< 5.0 %
Sulphated ash	< 15.0 %
Suitability for microbiolog	Passes test
Heat resistant bacteria	Not detected