

Technical Datasheet

Article No. 9257

**Yeast Extract for microbiology
(Powder)**

For Laboratory use.

Parameter	Specification
Amino Nitrogen (AN)	min. 4.5 %
Total Nitrogen (TN)	min. 10 %
Loss on drying	max. 5.0 %
Sodium Chloride (NaCl)	max. 0.5 %
pH (2 %)	6.0–7.2
Standard Plate Count	max. 5,000 col/g
Molds and yeast	max. 50 col/g
Coliforms	negative
Salmonella	negative