

Technical Datasheet

Article No. DG1250

Formic acid puriss.
98–100% HCOOH

For laboratory use.

Parameter	Value
Appearance	Clear, colourless liquid
Boiling point	101 °C
Density (20 °C)	1.220–1.221 g/cm ³
Mol.-Weight	46.026 g/mol
Assay (titr.)	min. 99.0 %
Acetic acid (CH ₃ COOH)	max. 0.005 %
Chloride (Cl)	max. 0.0001 %
Sulfate (SO ₄)	max. 0.0001 %
Iron (Fe)	max. 0.0001 %
Other heavy metals	max. 0.0001 %
Residue on evaporation	max. 0,002 %